



“The Gold Standard for Dicom Image Viewing and Diagnosis”



eFilm™

System Requirements

Note: A faster processor and more powerful video cards will improve performance.

Supported Operating Systems

- Windows 7-Home Premium, Professional, Enterprise and Ultimate (32 or 64-bit)
- Windows Vista - Enterprise, Business and Home Premium Editions
- Windows XP Professional (SP2)

Required Hardware

To guarantee optimum performance of eFilm™, your computer system should meet the specifications given below.

To get the most out of eFilm™, Merge Healthcare recommends that your computer meet the following specifications:

- Pentium 4-class processor
- 4 GB RAM
- 40 GB free hard drive space (1 GB for installing eFilm and 39 GB for image storage)
- Display resolution 1280 x 1024 with 32 bit/true color

For 3D Volume Rendering

- Video Card with at least 128 MB video
- RAM that fully supports DirectX 8.1 or higher

YOUR IMAGE IS OUR BUSINESS.



a VCA ANTECH company

800.268.5354 soundeklin.com